
BML-284

产品编号: D60909

CAS: 853220-52-7

分子式: C₁₉H₁₈N₄O₃

纯度: 98%

InChi Key: InChIKey=FABQUVYDAXWUQP-UHFFFAOYSA-N

Smiles: COc1cccc(c1)-c1cc(NCc2ccc3OCOc3c2)nc(N)n1

外观: 固体粉末

作用通路: Cytoskeletal Signaling

作用靶点: Wnt/beta-catenin

溶解性: Soluble in DMSO

保存条件: -20°C

产品介绍: BML-284 is a potent and cell-permeable Wnt signaling activator. BML-284 induces TCF-dependent transcriptional activity with an EC₅₀ of 700 nM