
TW-37

产品编号: D60923

CAS: 877877-35-5

分子式: C₃₃H₃₅NO₆S

纯度: 96%

InChi Key: InChIKey=PQAPVTKIEGUPRN-UHFFFAOYSA-N

Smiles: CC(C)c1cccc1Cc1c(c(c(c1)C(=O)Nc1ccc(cc1)S(=O)(=O)c1cccc1C(C)(C)O)O)O

外观: 固体粉末

作用通路: Apoptosis

作用靶点: BCL

溶解性: DMSO: ≥ 42 mg/mL (73.21 mM)

保存条件: -20°C

产品介绍: TW-37 is a potent Bcl-2 inhibitor with Ki values of 260, 290 and 1110 nM for Mcl-1, Bcl-2 and Bcl-xL, respectively.