
HJC0350

产品编号: D60925

CAS: 885434-70-8

分子式: C₁₅H₁₉NO₂S

纯度: 98%

InChi Key: InChIKey=AFZWZVLPIMHLSE-UHFFFAOYSA-N

Smiles: Cc1cc(c(c(c1)C)S(=O)(=O)n1cc(cc1C)C)C

外观: 固体粉末

作用通路: GPCR/GProtein

作用靶点: cAMP

溶解性: DMSO: 33.33 mg/mL (120.16 mM; Need ultrasonic)

保存条件: RT

产品介绍: HJC0350 is a potent and specific EPAC2 antagonist with an IC₅₀ of 0.3 μM.