
PF-562271 besylate

产品编号: D60954

CAS: 939791-38-5

分子式: C₂₁H₂₀F₃N₇O₃S · C₆H₆O₃S

纯度: 98%

InChi Key: InChIKey=LKLWTLXTOVZF AE-UHFFFAOYSA-N

Smiles: OS(=O)(=O)c1ccccc1.CN(c1c(CNc2c(cnc(Nc3ccc4NC(=O)Cc4c3)n2)C(F)(F)F)cccn1)S(=O)(=O)C

外观: 固体粉末

作用通路: Angiogenesis

作用靶点: CDK; FAK; PYK2

溶解性: Soluble in DMSO

保存条件: -20°C

产品介绍: PF-562271 (VS-6062) besylate is a potent ATP-competitive, reversible inhibitor of FAK and Pyk2 kinase, with an IC₅₀ of 1.5 nM and 13 nM, respectively