
JNJ-38877605

产品编号: D60959

CAS: 943540-75-8

分子式: C₁₉H₁₃F₂N₇

纯度: 99%

InChi Key: InChIKey=JRWCBEOAFGHNNU-UHFFFAOYSA-N

Smiles: Cn1cc(cn1)-c1ccc2nnc(n2n1)C(F)(F)c1ccc2ncccc2c1

外观: 固体粉末

作用通路: Tyrosine Kinase/Adaptors

作用靶点: c-Met/HGFR

溶解性: H₂O:<1 mg/mL Ethanol:<1 mg/mL DMSO:35 mg/mL (92.8 mM)

保存条件: -20°C

产品介绍: JNJ-38877605 is an ATP-competitive c-Met inhibitor (IC₅₀: 4 nM), 600-fold selective for c-Met than 200 other tyrosine and serine-threonine kinases.