
AMG 900

产品编号: D60962

CAS: 945595-80-2

分子式: C₂₈H₂₁N₇O₅

纯度: 99%

InChi Key: InChIKey=IVUGFMLRJOCGAS-UHFFFAOYSA-N

Smiles: CC1=CSC(=C1)C1=NN=C(NC2=CC=C(OC3=C(C=CC=N3)C3=NC(N)=NC=C3)C=C2)C2=C1C=CC=C2

外观: 固体粉末

作用通路: Cell Cycle/Checkpoint

作用靶点: Aurora Kinase; p38 MAPK; Tyrosine Kinases

溶解性: DMSO: ≥ 50 mg/mL (99.29 mM)

保存条件: 2-8°C

产品介绍: AMG 900 is a potent and highly selective pan-Aurora kinases inhibitor with IC₅₀ of 5 nM, 4 nM and 1 nM for Aurora A, B and C, respectively.