
PI-3065

产品编号: D60966

CAS: 955977-50-1

分子式: C₂₇H₃₁N₆O₃

纯度: 98%

InChi Key: InChIKey=YDNOHCOYQVZOMC-UHFFFAOYSA-N

Smiles: c1(c2c(nc(n1)c1c(ccc3c1cc[nH]3)F)cc(s2)CN1CCN(CC1)CC1CC1)N1CCOCC1

外观: 固体粉末

作用通路: DNA Damage/DNA Repair

作用靶点: DNA-PK; mTOR; PI3K

溶解性: Soluble in DMSO

保存条件: -20°C

产品介绍: PI-3065 is a novel potent and selective PI3K p110δ inhibitor.