
BI-847325

产品编号: D60998

CAS: 1207293-36-4

分子式: C₂₉H₂₈N₄O₂

纯度: 95%

InChi Key: InChIKey=FLBNLJLONKAPLR-DQSJHHFOSA-N

Smiles: CCNC(=O)C#Cc1ccc2c(NC(=O)\C2=C(/Nc2ccc(CN(C)C)cc2)c2ccccc2)c1

外观: 固体粉末

作用通路: Apoptosis

作用靶点: Apoptosis; Aurora Kinase; MEK

溶解性: Soluble in DMSO、Ethanol

保存条件: -20°C

产品介绍: BI-847325 is a selective dual inhibitor of MEK and aurora kinases (AK) with IC₅₀ values of 4 and 15 nM for human MEK2 and AK-C, respectively.