

---

## TAK-632

产品编号: D61004

CAS: 1228591-30-7

分子式: C<sub>27</sub>H<sub>18</sub>F<sub>4</sub>N<sub>4</sub>O<sub>3</sub>S

纯度: 98%

InChi Key: InChIKey=OJFKUJDRGJSAQB-UHFFFAOYSA-N

Smiles: FC1=C(NC(=O)CC2=CC(=CC=C2)C(F)(F)F)C=C(OC2=C(C#N)C3=C(C=C2)N=C(NC(=O)C2CC2)S3)C=C1

外观: 固体粉末

作用通路: Angiogenesis

作用靶点: Aurora Kinase; FGFR; PDGFR; Raf

溶解性: DMSO: 100 mg/mL (180.34 mM; Need ultrasonic)

保存条件: 2-8°C

产品介绍: TAK-632 is a potent pan-RAF inhibitor with IC<sub>50</sub> of 1.4, 2.4 and 8.3 nM for CRAF, BRAFV600E, BRAFWT, respectively.