
MG 149

产品编号: D61010

CAS: 1243583-85-8

分子式: C₂₂H₂₈O₃

纯度: 97%

InChi Key: InChIKey=WBHQYBZRTAEHRR-UHFFFAOYSA-N

Smiles: Oc1c(C(=O)O)c(CCc2ccc(CCCCCC)cc2)ccc1

外观: 固体粉末

作用通路: Apoptosis

作用靶点: Apoptosis; Epigenetic Reader Domain; Histone Acetyltransferase

溶解性: Soluble in DMSO

保存条件: -20°C

产品介绍: MG149 (Tip60 HAT inhibitor) is a selective and potent Tip60 inhibitor with IC₅₀ of 74 uM, similar potency for MOF (IC₅₀= 47 uM); little potent for PCAF and p300 (IC₅₀ >200 uM)