

---

## CP21R7

产品编号: D61011

CAS: 125314-13-8

分子式: C<sub>19</sub>H<sub>15</sub>N<sub>3</sub>O<sub>2</sub>

纯度: 98%

InChi Key: InChIKey=RGTAEYDIDMGJLX-UHFFFAOYSA-N

Smiles: C1(=O)C(=C(C(=O)N1)c1cc(ccc1)N)c1cn(c2c1cccc2)C

外观: 固体粉末

作用通路: Chromatin/Epigenetic

作用靶点: GSK-3; PKC

溶解性: Soluble in DMSO

保存条件: -20°C

产品介绍: CP21R7 is potent GSK-3 $\beta$  inhibitor, with an IC<sub>50</sub> of 1.8 nM; CP21R7 also shows inhibitory activity against PKC $\alpha$ , with an IC<sub>50</sub> of 1900 nM.