
EPZ004777

产品编号: D61022

CAS: 1338466-77-5

分子式: C₂₈H₄₁N₇O₄

纯度: 99%

InChi Key: InChIKey=WXRGFPHDRFQODR-ICLZECGLSA-N

Smiles: CC(C)N(CCCNC(=O)Nc1ccc(cc1)C(C)(C)C[C@H]1O[C@H]([C@H](O)[C@@H]1O)n1ccc2c(N)ncnc12

外观: 固体粉末

作用通路: Apoptosis

作用靶点: Apoptosis; Histone Methyltransferase

溶解性: Soluble in DMSO、Ethanol

保存条件: -20°C

产品介绍: EPZ004777 is a potent and selective DOT1L inhibitor with an IC₅₀ of 0.4 nM.