

---

## GSK503

产品编号: D61023

CAS: 1346572-63-1

分子式: C<sub>31</sub>H<sub>38</sub>N<sub>6</sub>O<sub>2</sub>

纯度: 99%

InChi Key: InChiKey=HRDQQHUKUIKFHT-UHFFFAOYSA-N

Smiles: CC(C)n1cc(C)c2c(cc(cc12)c1cnc(cc1)N1CCN(C)CC1)C(=O)NCc1c(C)cc(C)[nH]c1=O

外观: 固体粉末

作用通路: Chromatin/Epigenetic

作用靶点: Histone Methyltransferase

溶解性: Soluble in DMSO

保存条件: -20°C

产品介绍: GSK503 is a potent and specific inhibitor of EZH2 methyltransferase with Kiapp values of 3 to 27 nM.