
Senexin A

产品编号: D61025

CAS: 1366002-50-7

分子式: C₁₇H₁₄N₄

纯度: 99%

InChi Key: InChIKey=XBJCNHGQFJFCOY-UHFFFAOYSA-N

Smiles: N#CC=1C=CC2=NC=NC(NCCC=3C=CC=CC3)=C2C1

外观: 固体粉末

作用通路: Cell Cycle/Checkpoint

作用靶点: CDK

溶解性: Soluble in DMSO、Ethanol

保存条件: 2-8°C

产品介绍: Senexin A is an effective and selective CDK8 inhibitor that also inhibits CDK19 with K_d values of 0.83 microns and 0.31 microns, respectively.