
MI-503

产品编号: D61065

CAS: 1857417-13-0

分子式: C₂₈H₂₇F₃N₈S

纯度: 98%

InChi Key: InChIKey=DETOMBLLEOZTMZ-UHFFFAOYSA-N

Smiles: Cc1c(CN2CCC(CC2)Nc2ncnc3sc(CC(F)(F)F)cc23)ccc2n(Cc3cn[nH]c3)c(cc12)C#N

外观: 固体粉末

作用通路: Chromatin/Epigenetic

作用靶点: Epigenetic Reader Domain; Histone Methyltransferase

溶解性: Soluble in DMSO

保存条件: -20°C

产品介绍: MI-503 is a highly effective and orally bioavailable inhibitor of the menin-mLL interaction.