
TCS JNK 5a

产品编号: D61082

CAS: 312917-14-9

分子式: C₂₀H₁₆N₂O₂S

纯度: 98%

InChi Key: InChIKey=WQGDQGAFSDMBLA-UHFFFAOYSA-N

Smiles: C1CCc2c(C1)c(c(s2)NC(=O)c1cccc2ccccc12)C#N

外观: 固体粉末

作用靶点: Apoptosis; JNK

溶解性: DMSO: 20 mg/mL (60.16 mM; Need ultrasonic)

保存条件: RT

产品介绍: TCS JNK 5a is a potent JNK3 inhibitor with a pIC₅₀ of 6.7. TCS JNK 5a also inhibits JNK2 with a pIC₅₀ of 6.5.