
Eprodinate disodium

产品编号: D61086

CAS: 36589-58-9

分子式: C₃H₆Na₂O₆S₂

纯度: 98%

InChi Key: InChIKey=DKGJFKPIUSHDIT-UHFFFAOYSA-L

Smiles: [Na+].[O-]S(=O)(=O)CCCS(=O)(=O)[O-].[Na+]

外观: 固体粉末

作用靶点: Others

溶解性: H₂O: 50 mg/mL (201.46 mM; Need ultrasonic) DMSO: 7.69 mg/mL (30.98 mM; Need ultrasonic)

保存条件: RT

产品介绍: Eprodinate disodium (NC-503) is the orally available disodium salt form of Eprodinate, a negatively charged sulfonated inhibitor of fibrillogenesis, that can be used in the treatment of amyloid A (AA) amyloidosis.