
WH-4-023

产品编号: D61115

CAS: 837422-57-8

分子式: C₃₂H₃₆N₆O₄

纯度: 99%

InChi Key: InChIKey=NBTNHSGBRGTFJS-UHFFFAOYSA-N

Smiles: COc1ccc(N(C(=O)Oc2c(C)cccc2C)c2nc(Nc3ccc(N4CCN(C)CC4)cc3)ncc2)c(OC)c1

外观: 固体粉末

作用通路: Angiogenesis

作用靶点: Src

溶解性: Soluble in DMSO

保存条件: -20°C

产品介绍: WH-4-023 is a potent and selective dual Lck/Src inhibitor with IC₅₀ of 2 nM/6 nM for Lck and Src kinase respectively; little inhibition on p38α and KDR.