

---

## PF-04217903

产品编号: D61132

CAS: 956905-27-4

分子式: C<sub>19</sub>H<sub>16</sub>N<sub>8</sub>O

纯度: ≥98%

InChi Key: InChIKey=PDMUGYOXRHVNMU-UHFFFAOYSA-N

Smiles: OCCn1cc(cn1)c1cnc2nnn(Cc3cc4c(cc3)nccc4)c2n1

外观: 固体粉末

作用通路: Tyrosine Kinase/Adaptors

作用靶点: c-Met/HGFR

溶解性: Soluble in DMSO

保存条件: -20°C

产品介绍: PF-04217903 is a potent ATP-competitive c-Met kinase inhibitor with K<sub>i</sub> of 4.8 nM for human c-Met. PF-04217903 shows more than 1,000-fold selectivity relative to 208 kinases.  
Antiangiogenic properties