
Tegobuvir

产品编号: D61136

CAS: 1000787-75-6

分子式: C₂₅H₁₄F₇N₅

纯度: 97%

InChi Key: InChIKey=XBEQSQDCBSKCHJ-UHFFFAOYSA-N

Smiles: Fc1ccccc1-c1nc2ccn(Cc3ccc(nn3)-c3ccc(cc3C(F)(F)F)C(F)(F)F)cc2n1

外观: 固体粉末

作用通路: Microbiology/Virology

作用靶点: HCV Protease

溶解性: DMSO: ≥ 50 mg/mL (96.64 mM)

保存条件: -20°C

产品介绍: Tegobuvir is a specific, covalent inhibitor of the HCV NS5B polymerase.