
GNE-477

产品编号: D61152

CAS: 1032754-81-6

分子式: C₂₁H₂₈N₈O₃S₂

纯度: 98%

InChi Key: InChIKey=AKKCGLXULFRAET-UHFFFAOYSA-N

Smiles: Cc1c(sc2c1nc(nc2N1CCOCC1)c1cnc(nc1)N)CN1CCN(CC1)S(=O)(=O)C

外观: 固体粉末

作用通路: PI3K/Akt/mTOR signaling

作用靶点: mTOR; PI3K

溶解性: Soluble in DMSO

保存条件: -20°C

产品介绍: GNE-477 is a potent and efficacious dual PI3K (IC₅₀=4 nM)/mTOR(K_i=21 nM) inhibitor.