
ST271

产品编号: D61156

CAS: 106392-48-7

分子式: C₁₆H₂₀N₂O₂

纯度: 98%

InChi Key: InChIKey=HMLQHJRJKQDTKW-LFYBBSHMSA-N

Smiles: CC(C)c1cc(cc(c1O)C(C)C)/C=C(\C#N)/C(=O)N

外观: 固体粉末

作用通路: Metabolism

作用靶点: Phospholipase; Tyrosine Kinases

溶解性: DMSO: ≥ 110 mg/mL (403.91 mM)

保存条件: 2-8°C

产品介绍: ST271 is a potent inhibitor of protein tyrosine kinase (PTK), inhibits phospholipase D activation stimulated by fMet-Leu-Phe and PAF, with IC₅₀s of 6.7 and 9 μM, respectively.