
Seratrodast

产品编号: D61165

CAS: 112665-43-7

分子式: C₂₂H₂₆O₄

纯度: 99%

InChi Key: InChIKey=ZBVKEHDGYSLCCC-UHFFFAOYSA-N

Smiles: CC1=C(C)C(=O)C(C(CCCCCC(O)=O)c2ccccc2)=C(C)C1=O

外观: 固体粉末

作用通路: DNA Damage/DNA Repair

作用靶点: COX; PPAR; Prostaglandin Receptor

溶解性: DMSO: ≥ 100 mg/mL (282.14 mM)

保存条件: -20°C

产品介绍: Seratrodast(AA 2414) is an anti-asthmatic agent and a potent and selective thromboxane A₂ receptor (TP) antagonist.