

---

## E-4031

产品编号: D61169

CAS: 113559-13-0

分子式: C<sub>21</sub>H<sub>29</sub>Cl<sub>2</sub>N<sub>3</sub>O<sub>3</sub>S

纯度: 98%

InChi Key: InChIKey=ZQBNWMFBOSOOOLX-UHFFFAOYSA-N

Smiles: Cl.Cl.Cc1cccc(CCN2CCC(CC2)C(=O)c2ccc(NS(C)(=O)=O)cc2)n1

外观: 固体粉末

作用通路: Angiogenesis

作用靶点: EGFR; Potassium Channel

溶解性: Soluble in H<sub>2</sub>O

保存条件: -20°C

产品介绍: E-4031 is a selective hERG potassium channel blocker for use in class III anti-arrhythmic studies