
SB290157 trifluoroacetate

产品编号: D61171

CAS: 1140525-25-2

分子式: C₂₄H₂₉F₃N₄O₆

纯度: ≥97%(HPLC)

InChi Key: InChIKey=ZJRMPPVJAQWGEG-FYZYNONXSA-N

Smiles: OC(=O)C(F)(F)F.NC(=N)NCCC[C@H](NC(=O)COCC(c1ccccc1)c1ccccc1)C(O)=O

外观: 固体粉末

作用通路: Immunology/Inflammation

作用靶点: Complement System

溶解性: Soluble in DMSO、Ethanol

保存条件: -20°C

产品介绍: SB290157 trifluoroacetate is a potent and selective antagonist of C3a receptor(IC₅₀ of 200 nM).