
EDO-S101

产品编号: D61225

CAS: 1236199-60-2

分子式: C₁₉H₂₈Cl₂N₄O₂

纯度: 98%

InChi Key: InChIKey=GISXTRIGVCKQBX-UHFFFAOYSA-N

Smiles: O=C(NO)CCCCCc1nc2cc(N(CCCl)CCCl)ccc2n1C

外观: 固体粉末

作用通路: Chromatin/Epigenetic

作用靶点: HDAC

溶解性: Soluble in DMSO

保存条件: -20°C

产品介绍: Tinostamustine (EDO-S101) is a pan HDAC inhibitor; inhibits HDAC6, HDAC1, HDAC2 and HDAC3 with IC₅₀ values of 6 nM, 9 nM, 9 nM and 25 nM, respectively