
CNX-1351

产品编号: D61243

CAS: 1276105-89-5

分子式: C₃₀H₃₅N₇O₃S

纯度: 99%

InChi Key: InChIKey=DLNUPKDFXMWRFP-UHFFFAOYSA-N

Smiles: C1CN(CCN1C(=O)CCC(=O)C=C(C)C)Cc1cc2c(s1)c(nc(n2)c1cccc2c1cn[nH]2)N1CCOCC1

外观: 固体粉末

作用通路: PI3K/Akt/mTOR signaling

作用靶点: PI3K

溶解性: Soluble in DMSO

保存条件: -20°C

产品介绍: CNX-1351 is a potent and isoform-selective targeted covalent PI3K α inhibitor with IC₅₀ of 6.8 nM.