
SR1001

产品编号: D61272

CAS: 1335106-03-0

分子式: C₁₅H₁₃F₆N₃O₄S₂

纯度: 99%

InChi Key: InChIKey=OZBSSKGBKHOLGA-UHFFFAOYSA-N

Smiles: CC(=O)Nc1nc(C)c(s1)S(=O)(=O)Nc1ccc(cc1)C(O)(C(F)(F)F)C(F)(F)F

外观: 固体粉末

作用通路: Metabolism

作用靶点: ROR

溶解性: Soluble in DMSO、Ethanol

保存条件: -20°C

产品介绍: SR 1001 is an inverse agonist for ROR α /ROR γ (K_i: 172/111 nM).