
ML240

产品编号: D61277

CAS: 1346527-98-7

分子式: C₂₃H₂₀N₆O

纯度: 99%

InChi Key: InChIKey=NHAMBLRUUJAFOY-UHFFFAOYSA-N

Smiles: COc1cccc2c1nc(nc2NCc1cccc1)n1c2cccc2nc1N

外观: 固体粉末

作用靶点: p97

溶解性: Soluble in DMSO

保存条件: -20°C

产品介绍: ML240 is a potent p97 inhibitor, inhibiting p97 ATPase with IC₅₀ value of 100 nM.