
Leniolisib

产品编号: D61282

CAS: 1354690-24-6

分子式: C₂₁H₂₅F₃N₆O₂

纯度: 99%

InChi Key: InChIKey=MWKYMZXCGYXLPL-ZDUSSCGKSA-N

Smiles: CCC(=O)N1CC[C@@H](C1)Nc1ncnc2CCN(Cc12)c1cnc(OC)c(c1)C(F)(F)F

外观: 固体粉末

作用通路: DNA Damage/DNA Repair

作用靶点: DNA-PK; PI3K

溶解性: Soluble in DMSO

保存条件: -20°C

产品介绍: Leniolisib (CDZ173) is a potent and selective PI3Kδ inhibitor. Leniolisib has the potential for immunodeficiency disorders treatment.