
SR-3029

产品编号: D61344

CAS: 1454585-06-8

分子式: C₂₃H₁₉F₃N₈O

纯度: 98%

InChi Key: InChIKey=CEBMEQPREMCWOL-UHFFFAOYSA-N

Smiles: C1COCCN1c1nc2c(c(n1)NCc1nc3c([nH]1)ccc(c3F)F)ncn2c1cc(ccc1)F

外观: 固体粉末

作用通路: Metabolism

作用靶点: Casein Kinase

溶解性: Soluble in DMSO

保存条件: -20°C

产品介绍: SR-3029 is a potent and ATP competitive CK1δ and CK1ε inhibitor, with IC₅₀s of 44 nM and 260 nM, respectively, and Kis of 97 nM for both kinases.