
ML355

产品编号: D61359

CAS: 1532593-30-8

分子式: C₂₁H₁₉N₃O₄S₂

纯度: 98%

InChi Key: InChIKey=OWHBVKBNRYMIN-UHFFFAOYSA-N

Smiles: COC1=CC=CC(=C1O)CNC2=CC=C(C=C2)S(=O)(=O)NC3=NC4=CC=CC=C4S3

外观: 固体粉末

作用通路: Metabolism

作用靶点: Lipoxygenase

溶解性: Soluble in DMSO

保存条件: -20°C

产品介绍: ML355, a specific, effective 12-Lipoxygenase(12-LOX) inhibitors(IC₅₀=0.34 μM), possess favorable ADME properties.