
ADU-S100 ammonium salt

产品编号: D61384

CAS: 1638750-96-5

分子式: C₂₀H₃₀N₁₂O₁₀P₂S₂

纯度: 99%

InChi Key: InChIKey=VZYXAGGQHVUTHM-ITQXCAEYSA-N

Smiles: NC1=C(C2=NC=N1)N=CN2[C@@H]3O[C@]([H])(CO4)[C@@]([H])(O[P@@](OC[C@@]5([C@]([H])([C@]([H])([C@@H](O5)N6C7=NC=NC(N)=C7N=C6)O[P@@]4(S)=O)O)[H])(S=O)[C@H]3O.N

外观: 固体粉末

作用通路: Immunology/Inflammation

作用靶点: STING

溶解性: Soluble in DMSO、Methanol、H₂O

保存条件: -20°C

产品介绍: ADU-S100 ammonium salt (MIW815 ammonium salt), an activator of stimulator of interferon genes (STING), leads to potent and systemic tumor regression and immunity