
Tenalisib

产品编号: D61385

CAS: 1639417-53-0

分子式: C₂₃H₁₈FN₅O₂

纯度: 98%

InChi Key: InChIKey=HDXDQPRPFRKGKZ-INIZCTEOSA-N

Smiles: CC[C@@H](C1=C(C(=O)C2=CC=CC=C2O1)C3=CC(=CC=C3)F)NC4=NC=NC5=C4NC=N5

外观: 固体粉末

作用通路: PI3K/Akt/mTOR signaling

作用靶点: PI3K

溶解性: Soluble in DMSO

保存条件: -20°C

产品介绍: Tenalisib (RP6530) is a novel, potent, and selective PI3Kδ and PI3Kγ inhibitor with IC₅₀ values of 25 and 33 nM, respectively.