
Saroglitazar Magnesium

产品编号: D61386

CAS: 1639792-20-3

分子式: C₂₅H₂₈NO₄S.1/2Mg

纯度: 98%

InChi Key: InChIKey=UJYFZCVPOSZDMK-YPPDDXJESA-L

Smiles: CCO[C@H](C([O-])=O)CC1=CC=C(OCCN2C(C)=CC=C2C3=CC=C(SC)C=C3)C=C1.CCO[C@H](C([O-])=O)CC4=CC=C(OCCN5C(C)=CC=C5C6=CC=C(SC)C=C6)C=C4.[Mg+2]

外观: 固体粉末

作用靶点: Others

溶解性: Soluble in DMSO

保存条件: -20°C

产品介绍: Saroglitazar magnesium is a novel peroxisome proliferator-activated receptor (PPAR) agonist with predominant PPAR α and moderate PPAR γ activity with EC₅₀ values of 0.65 pM and 3 nM in HepG2 cells, respectively.