
Mexenone

产品编号: D61387

CAS: 1641-17-4

分子式: C₁₅H₁₄O₃

纯度: 98%

InChi Key: InChIKey=MJVGBKJNTFCUJM-UHFFFAOYSA-N

Smiles: COc1cc(O)c(cc1)C(=O)c1ccc(C)cc1

外观: 固体粉末

作用靶点: Others

溶解性: Soluble in DMSO

保存条件: -20°C

产品介绍: Mexenone is a potent benzophenone-type UV filter. Mexenone is used for sunscreensing agent