
S63845

产品编号: D61417

CAS: 1799633-27-4

分子式: C₃₉H₃₇ClF₄N₆O₆S

纯度: 99%

InChi Key: InChIKey=ZFBHXVOCZBPADE-SSEXGKCCSA-N

Smiles: O=C(O)[C@@H](CC1=C(OCC2=CC=NN2CC(F)(F)F)C=CC=C1)OC3=C4C(SC(C5=CC=C(F)O5)=[C@@]4[C@@]6=C(C)C(Cl)=C(OCCN7CCN(C)CC7)C=C6)=NC=N3

外观: 固体粉末

作用通路: Apoptosis

作用靶点: Bcl

溶解性: Soluble in DMSO

保存条件: -20°C

产品介绍: S63845 is a potent and selective myeloid cell leukemia 1 (MCL1) inhibitor with a K_d of 0.19 nM for human MCL1