
PI4KIIIbeta-IN-10

产品编号: D61428

CAS: 1881233-39-1

分子式: C₂₂H₂₅N₃O₅S₂

纯度: 98%

InChi Key: InChiKey=PLUYFBRIGUAKBR-UHFFFAOYSA-N

Smiles: COC1=CC=C(C2=C(C)N=C(NC(C(C)(C)C)=O)S2)C=C1S(NC3=CC=C(O)C=C3)(=O)=O

外观: 固体粉末

作用通路: PI3K/Akt/mTOR signaling

作用靶点: PI3K; PI4K

溶解性: Soluble in DMSO

保存条件: -20°C

产品介绍: PI4KIIIbeta-IN-10 is a potent PI4KIIIβ inhibitor with an IC₅₀ of 3.6 nM.