

---

## Tipifarnib (S enantiomer)

产品编号: D61435

CAS: 192185-71-0

分子式: C<sub>27</sub>H<sub>22</sub>Cl<sub>2</sub>N<sub>4</sub>O

纯度: 99%

InChi Key: InChIKey=PLHJCIYEEKOWNM-MHZLTWQESA-N

Smiles: Cn1c([C@](c(cc2ccc2Cl)(c(cc2)cc(C(c3cccc(Cl)c3)=C3)c2N(C)C3=O)N)cnc1

外观: 固体粉末

作用通路: Metabolism

作用靶点: Transferase

溶解性: Soluble in DMSO

保存条件: -20°C

产品介绍: Tipifarnib S enantiomer is the S-enantiomer of Tipifarnib. Tipifarnib is a potent and specific farnesyltransferase (FTase) inhibitor with IC<sub>50</sub> of 0.6 nM. Tipifarnib S enantiomer is the less active isomer.