
WHI-P97

产品编号: D61451

CAS: 211555-05-4

分子式: C₁₆H₁₃Br₂N₃O₃

纯度: 99%

InChi Key: InChIKey=YVCXQRVWNQMZEI-UHFFFAOYSA-N

Smiles: Brc1cc(cc(Br)c1O)Nc3ncnc2cc(OC)c(cc23)OC

外观: 固体粉末

作用通路: Angiogenesis

作用靶点: JAK

溶解性: DMSO: 3.33 mg/mL (7.32 mM; ultrasonic and warming and heat to 80°C)

保存条件: -20°C

产品介绍: WHI-P97 is a potent and selective JAK-3 inhibitor. WHI-P97 is effective in preventing the development allergic asthma in vivo