
TCS 401

产品编号: D61475

CAS: 243966-09-8

分子式: C₁₀H₁₁ClN₂O₅S

纯度: 98%

InChi Key: InChIKey=LQGCAMWQDSYOAY-UHFFFAOYSA-N

Smiles: Cl.OC(=O)C(=O)Nc1sc2CNCCc2c1C(O)=O

外观: 固体粉末

作用通路: Metabolism

作用靶点: Phosphatase

溶解性: Soluble in DMSO

保存条件: -20°C

产品介绍: TCS 401 is a selective inhibitor of protein-tyrosine phosphatase 1B (PTP1B; K_i = 0.29 μM)