
Toltrazuril (sulfone)

产品编号: D61651

CAS: 69004-04-2

分子式: C₁₈H₁₄F₃N₃O₆S

纯度: 98%

InChi Key: InChIKey=VBUNOIXRZJNAD-UHFFFAOYSA-N

Smiles: Cc1cc(ccc1Oc1ccc(cc1)S(=O)(=O)C(F)(F)F)-n1c(=O)[nH]c(=O)n(C)c1=O

外观: 固体粉末

作用通路: Microbiology/Virology

作用靶点: Parasite

溶解性: Soluble in DMSO

保存条件: 2-8°C

产品介绍: Toltrazuril sulfone (Ponazuril) is a metabolite of Toltrazuril (HY-B0175), with antiprotozoal activity. Toltrazuril sulfone is a triazine anticoccidial that is developed to prevent coccidiosis in poultry