
Bipenquinatate

产品编号: D61829

CAS: 96201-88-6

分子式: C₂₃H₁₄F₂NNaO₂

纯度: 98%

InChi Key: InChIKey=PZOHOALJQOFNTB-UHFFFAOYSA-M

Smiles: [Na+].Cc1c(nc2ccc(F)cc2c1C([O-])=O)-c1ccc(cc1)-c1ccccc1F

外观: 固体粉末

作用通路: Metabolism

作用靶点: Dehydrogenase

溶解性: Soluble in DMSO

保存条件: -20°C

产品介绍: Brequinar (sodium) is a biochemical reagent that can be used as a biological material or organic compound for life science related research