
ICA 110381

产品编号: D61844

CAS: 325457-99-6

分子式: C₁₂H₈Cl₂N₂O

纯度: 98%

InChi Key: InChIKey=FOYOGIUAMYSVIT-UHFFFAOYSA-N

Smiles: O=C(Nc1ccc(Cl)nc1)c1ccc(Cl)cc1

外观: 固体粉末

作用通路: Membrane transporter/Ion channel

作用靶点: Potassium Channel

溶解性: Soluble in DMSO、Ethanol

保存条件: -20°C

产品介绍: ICA 110381 is a KCNQ2/Q3 potassium channel opener for the treatment of epilepsy