
ST034307

产品编号: D61969

CAS: 133406-29-8

分子式: C₁₀H₄Cl₄O₂

纯度: 98%

InChi Key: InChIKey=NTDHYMSVCBGQJF-UHFFFAOYSA-N

Smiles: ClC1=CC2=C(OC(=CC2=O)C(Cl)(Cl)Cl)C=C1

外观: 固体粉末

作用通路: Neuroscience

作用靶点: Adenylyl cyclase

溶解性: Soluble in DMSO

保存条件: -20°C

产品介绍: ST034307 is a potent and selective adenylyl cyclase 1 (AC1) inhibitor, with IC₅₀ of 2.3 μM.