
ZK824190

产品编号: D61972

CAS: 2254001-81-3

分子式: C₂₃H₂₂F₂N₂O₄

纯度: 97%

InChi Key: InChIKey=ZDHMLADHONAIES-LJQANCHMSA-N

Smiles: CC[C@H](C(O)=O)Oc(c(F)c1nc(Oc2cc(-c3cc(CN)ccc3)ccc2)c1F

外观: 固体粉末

作用通路: Proteases/Proteasome

作用靶点: Serine Protease

溶解性: Soluble in DMSO、Ethanol

保存条件: 2-8°C

产品介绍: ZK824190 is an orally available and selective inhibitor of urokinase plasminogen activator (uPA). The IC₅₀ values are 237, 1600 and 1850 nM for uPA, tPA, and Plasmin, respectively.