
DSP-2230

产品编号: D61973

CAS: 1233231-30-5

分子式: C₂₀H₂₀F₃N₅O₂

纯度: 98%

InChi Key: InChIKey=HHXCJIMPEJSJTG-SNVBAGLBSA-N

Smiles: C[C@H](C(N)=O)NCc(n1C2CCC2)nc(cc2)c1nc2Oc(cc1F)cc(F)c1F

外观: 固体粉末

作用通路: Membrane transporter/Ion channel

作用靶点: Sodium Channel

溶解性: Soluble in DMSO

保存条件: -20°C

产品介绍: DSP-2230 is a selective Nav1.7/Nav1.8 blocker