

---

## Linrodostat

产品编号: D61974

CAS: 1923833-60-6

分子式: C<sub>24</sub>H<sub>24</sub>ClFN<sub>2</sub>O

纯度: 98%

InChi Key: InChIKey=KRTIYQIPSAGSBP-ZACQAIPSSA-N

Smiles: C[C@H]([C@H]1CC[C@H](CC1)c1ccnc2ccc(F)cc12)C(=O)Nc1ccc(Cl)cc1

外观: 固体粉末

作用通路: Metabolism

作用靶点: IDO; Indoleamine 2,3-Dioxygenase(IDO)

溶解性: Soluble in DMSO

保存条件: -20°C

产品介绍: Linrodostat (BMS-986205) is a selective and irreversible indoleamine 2,3-dioxygenase 1 (IDO1) inhibitor with an IC<sub>50</sub> value of 1.1 nM in IDO1-HEK293 cells. Linrodostat is well tolerated with potent pharmacodynamic activity in advanced cancers