
ML 297

产品编号: D61979

CAS: 1443246-62-5

分子式: C₁₇H₁₄F₂N₄O

纯度: 98%

InChi Key: InChIKey=IEKSMUSSYJUQMY-UHFFFAOYSA-N

Smiles: Fc1ccc(cc1F)NC(=O)Nc3cc(C)nn3c2ccccc2

外观: 固体粉末

作用通路: Membrane transporter/Ion channel

作用靶点: Potassium Channel

溶解性: Soluble in DMSO

保存条件: -20°C

产品介绍: ML 297 (VU 0456810) is a potent and selective GIRK1/2 activator, with an EC₅₀ of 0.16 μM. ML 297 is potential for the treatment of epilepsy