
NSC 228155

产品编号: D61984

CAS: 113104-25-9

分子式: C₁₁H₆N₄O₄S

纯度: 98%

InChi Key: InChIKey=ICCFXXDUYSPKOL-UHFFFAOYSA-N

Smiles: [O-][N+](=O)c1ccc(Sc2cccc[n+]2[O-])c2nonc12

外观: 固体粉末

作用通路: Angiogenesis

作用靶点: DNA/RNA Synthesis; EGFR; Epigenetic Reader Domain; Histone Acetyltransferase

溶解性: Soluble in DMSO

保存条件: -20°C

产品介绍: NSC228155 is an activator of EGFR, binding to the sEGFR dimerization domain II and modulate EGFR tyrosine phosphorylation.